

EASY MAINTENANCE Operator Guide

MOCOM MILLENNIUM BP, BM, BG, BH, BL Chamber Filter Location Version 1


Regular Operator maintenance is **essential** for good performance. Refer to **Appendix C - Maintenance** in your User Manual for a details specific to your model.

Recommended general operator maintenance is listed below:

<p>Daily</p>	<p>Clean the ‘rubber’ chamber seal (1) and stainless steel chamber face (2) against which the rubber seal sits. Use a clean soft lint free cloth (not paper) with mild detergent and warm water. Clean any debris and rinse well to remove any detergent residue. Dry the surfaces.</p> <p>Only if stubborn build-up is present on the chamber face, use a nylon scourer (eg. green Scotchbrite pad). Rinse well to remove any residue. Dry the surfaces.</p>
<p>Weekly</p>	<p>Clean the chamber rack (5), inner chamber (3), door plate (4), and trays. Remove the trays and chamber rack from the chamber. Clean as above, rinse to remove any residue and dry all surfaces.</p> <p>Clean the outer surfaces. Clean, rinse, removing any residue and dry.</p> <p>Clean the mesh filter or replace if necessary. Open the door of the steriliser and remove the plastic filter cover plate (6) with your fingers. Using a large flat blade screwdriver, carefully remove the filter holder (7), then take the mesh filter (8) out of the holder (this is a tight fit).</p> <p>Clean the mesh filter under running water. If necessary use a pointed tool to remove fluff and lint from inside. If necessary use a wire brush to remove build-up on the outside.</p> <p>If the mesh filter cannot be cleaned acceptably, replace with a new one (2MOC-0043).</p> <p>Replace the mesh filter in the filter holder and screw gently back into position. Add a little silicone lubricant to the o-ring if needed. Replace plastic filter cover plate.</p>
<p>As Needed</p>	<p>The chamber seal, bacterial filter and printer paper need to be replaced periodically—the frequency depends on usage. These items can be ordered from your local Mocom service provider.</p> <p>2MOC-0011 Door Seal, grey 2MOC-0020 Bacterial Filter 2MOC-0042 Chamber Mesh Filter 9PAPER-03 Thermal Paper Roll (9 rolls)</p>

Key:

1. ‘Rubber’ chamber seal
2. Stainless steel chamber face
3. Inner chamber
4. Door plate
5. Chamber rack
6. Plastic filter cover plate
7. Filter holder
8. Mesh filter
9. Bacterial filter



diagrams copyright © MOCOM srl.